

## **CoAxial** Spray Kit

Simplifying the Application of  
Platelet-Rich Plasma onto Graft Material

**BIOMET**  
BIOLOGICS

## One Surgeon. One Patient.

**Over 1 million times per year, Biomet helps one surgeon provide personalized care to one patient.**

The science and art of medical care is to provide the right solution for each individual patient. This requires clinical mastery, a human connection between the surgeon and the patient, and the right tools for each situation.

At Biomet, we strive to view our work through the eyes of one surgeon and one patient. We treat every solution we provide as if it's meant for a family member.

Our approach to innovation creates real solutions that assist each surgeon in the delivery of durable personalized care to each patient, whether that solution requires a minimally invasive surgical technique, advanced biomaterials or a patient-matched implant.

**When one surgeon connects with one patient to provide personalized care, the promise of medicine is fulfilled.**

# CoAxial Spray Kit

---

The Biomet Biologics CoAxial spray applicator provides an efficient and convenient way to load and apply platelet-rich plasma (PRP).

## CoAxial Spray Kit Advantages

- Two applicator syringes and two tips included in each kit allow user to completely load both applicators for time savings (Figure 1).
- Replaces current method of transferring PRP into a plastic cup on the sterile field.
- Simplistic design allows for loading of PRP in a single step.
- Ensures 10:1 ratio of liquids (Figure 2).
- Provides fine mist spray without need for additional equipment or air source (Figure 3).



Figure 1  
Contents of CoAxial  
Spray Kit (800-0270)



Figure 2




Figure 3

# CoAxial Spray Kit

## Instructions for Use



### Step One

Transfer the CoAxial Spray Kit and Tips to the sterile field. Pre-assemble Y-connector to each CoAxial sprayer.

**Note:** Ensure Y-connector is fully tightened.

### Step Two

Remove red cap on Y-connector and attach PRP-filled syringe.

**Note:** Ensure PRP-filled syringe is attached to red, angled end of Y-connector.



### Step Three

Remove clear cap on Y-connector and attach syringe containing clotting factors.

**Note:** Ensure PRP syringe is loaded on red connector side and clotting factors syringe is attached on clear, straight side.


### Step Four

Push down on the PRP-filled syringe plunger. Once the PRP-filled syringe plunger is pushed, the clotting factors syringe automatically loads into a separate chamber. This ensures a 10:1 ratio.

**Note:** Filling is complete when either the smaller or larger syringe is void of liquid. At this time, stop depressing plunger to ensure a 10:1 ratio.

# CoAxial Spray Kit

## Instructions for Use



### Step Five

Remove the Y-connector and attached syringes as one complete fixture from the CoAxial Sprayer.



### Step Six

Assemble the sterile spray tip onto the CoAxial Sprayer by twisting onto the threaded hub.

**Note:** If needed, the second CoAxial Sprayer can be loaded in the same fashion.



## Step Seven

The PRP can now be sprayed onto the graft material of choice.\*

\* The CoAxial Applicator System is designed to facilitate pre-mixing of allograft, autograft, or synthetic bone graft materials for application to an orthopedic surgical site, with autologous blood, plasma, platelet-rich plasma, or other specific blood component(s) as deemed necessary by the clinical use requirements. In addition, the CoAxial Applicator System is intended for the application of autologous blood, plasma, platelet-rich plasma, or other specific blood component(s), as deemed necessary by the surgeon's determination of the clinical use requirements to facilitate the preparation of soft tissue autograft or allograft material prior to the application of the graft material to a repair site.

## Ordering Information

Part Number	Description
800-0270	CoAxial Spray Kit with two spray tips
800-0274	CoAxial 4-inch Malleable Dual Spray Tip

All content herein is protected by copyright, trademarks and other intellectual property rights owned by or licensed to Biomet Inc. or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Biomet.

This material is intended for health care professionals and the Biomet sales force. Distribution to any other recipient is prohibited.

For complete product information, including indications, contraindications, warnings, precautions, and potential adverse effects, see the package insert and [www.Biomet.com](http://www.Biomet.com).



**Responsible Manufacturer**

Biomet Biologics, Inc.  
P.O. Box 587  
56 E. Bell Drive  
Warsaw, Indiana 46581-0587  
USA