

# The Rationale Behind the I.D.E.A.L.<sup>TM</sup> Femoral Tunnel Position Philosophy

by Stephen Howell, M.D.

A general consensus from *in vitro* and *in vivo* studies is that the ideal position of the femoral tunnel is when the ACL graft is **I**sometric, in the **D**irect fibers, **E**quidistant and **E**ccentric, **A**natomic, and **L**ow in tension. The **I.D.E.A.L.** acronym is the foundation for locating the position of the femoral guide wire intraoperatively in the transtibial, anteromedial portal, and two-incision surgical techniques.

The intraoperative check for correct positioning of the femoral tunnel is placement of the guide wire equidistant from the top to the bottom of the notch within the green zone and far enough posterior to position the femoral tunnel with a  $\leq 2$  mm tunnel backwall' (Figure 1 & 2). Centering the femoral guide wire in the green zone has the advantages of setting graft length isometrically,<sup>2</sup> anatomically centering the graft in the direct fibers of the origin of the native ACL,<sup>3</sup> and restoring low graft tension;<sup>2,4,5</sup> which are associated with high function, full motion, and stability.<sup>1</sup>

In contrast, positioning the femoral guide wire so the boundary of the femoral tunnel extends outside the green zone has three undesirable consequences:

1. Placement 'vertical' or above the green zone may cause high graft tension from impingement of the graft against the PCL in flexion, loss of flexion, and a high failure rate.<sup>5</sup>
2. Placement 'low' or below the green zone may cause high graft tension from non-isometric placement in extension, loss of extension, and a high failure rate.<sup>4,6</sup>
3. Placement anterior, so the femoral tunnel has  $>2$  mm backwall, may cause high graft tension from non-isometric placement in flexion, and loss of flexion.<sup>2</sup>

The advantage of the green zone with the I.D.E.A.L. philosophy is it allows an anterior medial or transtibial drilling technique, along with the ability to consistently place the tunnel within the green zone and allows a small amount of latitude in order to support the individuality of the patient's anatomy and notch.



Femoral Guide Wire



Femoral Tunnel



$< 2$  mm of Backwall

Figure 1

Figure 2

● Green Zone ● Tunnel Too Vertical or Low

## References

1. Smith, C.K.; Howell, S.M.; and Hull, M.L.: Anterior laxity, slippage, and recovery of function in the first year after tibialis allograft anterior cruciate ligament reconstruction. *Am J Sports Med*, 39(1): 78–88, 2011.
2. Zavras, T.D.; Race, A.; and Amis, A.A.: The effect of femoral attachment location on anterior cruciate ligament reconstruction: graft tension patterns and restoration of normal anterior-posterior laxity patterns. *Knee Surg Sports Traumatol Arthrosc*, 13(2): 92–100, 2005.
3. Sasaki, N., M.D.; *et al.*: The femoral insertion of the anterior cruciate ligament: Discrepancy between macroscopic and histological observations. *Arthroscopy*, 28(8): 1135–1146, 2012.
4. Markolf, K.L.; Park, S.; Jackson, S.R.; and McAllister, D.R.: Anterior-posterior and rotatory stability of single and double-bundle anterior cruciate ligament reconstructions. *J Bone Joint Surg Am*, 91(1): 107–18, 2009.
5. Simmons, R.; Howell, S.M.; and Hull, M.L.: Effect of the angle of the femoral and tibial tunnels in the coronal plane and incremental excision of the posterior cruciate ligament on tension of an anterior cruciate ligament graft; an in vitro study. *J Bone Joint Surg Am*, 85-A(6): 1018–29, 2003.
6. Rahr-Wagner, L.; Thillemann, T.M.; Pedersen, A.B; and Lind, M.C.: Increased risk of revision after anteroedial compared with transtibial drilling of the femoral tunnel during primary anterior cruciate ligament reconstruction: results from the Danish Knee Ligament Reconstruction Register. *Arthroscopy*, 29(1): 98–105, 2013.

This positioning philosophy was prepared by a licensed health care professional. Zimmer Biomet does not practice medicine and does not recommend any particular positioning for use on a specific patient. The surgeon is responsible for determining the appropriate positioning for each individual patient.

This material is intended for health care professionals and the Zimmer Biomet sales force only. Distribution to any other recipient is prohibited. All content herein is protected by copyright, trademarks and other intellectual property rights owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated. This material must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

Check for country product clearances and reference product specific instructions for use. For complete product information, including indications, contraindications, warnings, precautions, and potential adverse effects, see the package insert and Zimmer Biomet's website.

Not for distribution in France

©2016 Zimmer Biomet



**ZIMMER BIOMET**  
Your progress. Our promise.™