



Cable-Ready® Cable Grip System Cable Buttons

Data Sheet



47-2232-060-01

47-2232-060-55

47-2232-060-35

Enhance intraoperative flexibility and surgeon choice

Cable-Ready Cable Grip Buttons

The *Cable-Ready*® Cable Grip System Button implants are designed to be used in conjunction with *NCB*® Polyaxial Locking, *NCB* Periprosthetic and *Zimmer*® Periarticular Locking Plates.

The Cable Button is threaded into a vacant screw hole of a Zimmer Locking Plate to provide a positioning point for the Cerclage Cable. The *NCB* Cable Button for *NCB* Polyaxial Locking Plate is made of Ti-6Al-4V Alloy and the *Zimmer* Periarticular Locking Plate Cable Buttons are made of 22-13-5 Stainless Steel. The buttons have a 2.5mm hex feature that mates with the standard hex screw driver.

Indications

The Cable Button is intended for use in combination with the Zimmer Locking Bone Plating Systems and Cerclage Cables, to stabilize multiple fractures or butterfly fragments in long bones.

Contraindications

- Infection
- Sternal or spinal fractures or defects

The following three buttons are available:

- Cable Buttons for *NCB* Polyaxial Locking Plates and all *NCB* Periprosthetic Plates, Titanium alloy, 2.5mm Hex Drive (REF 47-2232-060-01)
- Cable Buttons for the 3.5mm *Zimmer* Periarticular Locking Humeral Plate, 3.5mm Proximal Lateral Tibial Plate, Distal Lateral Tibial Plate, Distal Medial Tibial Plate, and the 3.5mm *Zimmer* Universal Locking Plates, Stainless Steel, 2.5mm Hex Drive (REF 47-2232-060-35)
- Cable Buttons for the 5.5mm *Zimmer* Periarticular Locking Femoral Plate and 5.5mm Proximal Lateral Tibial Plate, Stainless Steel, 2.5mm Hex Drive (REF 47-2232-60-55)

Cable-Ready Cable Buttons

Prod. No.	Description	Size
47-2232-060-01	<i>NCB</i> Cable Button For <i>NCB</i> Polyaxial Locking Plate, Single, Titanium Alloy	2.5mm Hex
47-2232-060-35	3.5mm Peri Locking Cable Button, Single, Stainless Steel	2.5mm Hex
47-2232-060-55	5.5mm Peri Locking Cable Button, Single, Stainless Steel	2.5mm Hex

Cable-Ready Cable Buttons Instructions

Select a Stainless Steel cable for use with Stainless Steel locking plates or a Cobalt Chrome cable for use with Titanium alloy locking plates.

1 Cable Button Insertion

Select the appropriate cable button, using the 2.5mm hex head screwdriver thread it into the plate hole. Do not fully tighten to allow the slots in the button to align with the cable.

2 Pass the Cable

The cable is inserted into the tip of the Cable Passer and threaded until the cable exits the shaft of the Cable Passer.

3 Cable Insertion

Pass the cable through the button and feed it through the crimp.

4 Tensioning the Cable

Tension the cable. After the cable is tightened, lock down the Retensioning Bit to temporarily hold the cable. If necessary, move the Tensioner to the next cable and tension it. Continue this procedure until all the cables you would like to use are tensioned. (For additional information refer to the cable tensioner surgical technique)

5 Crimping the Cable

Carefully unlock the Retensioning Bit and re-tension the first cable. Crimp the first cable with the screwdriver.

6 Remove the Excess Cable

After the set screw is firmly seated, cut the excess cable flush with the connector body.



Disclaimer

This documentation is intended exclusively for physicians and is not intended for laypersons. Information on the products and procedures contained in this document is of a general nature and does not represent and does not constitute medical advice or recommendations. Because this information does not purport to constitute any diagnostic or therapeutic statement with regard to any individual medical case, each patient must be examined and advised individually, and this document does not replace the need for such examination and/or advice in whole or in part. Please refer to the package inserts for important product information, including, but not limited to, contraindications, warnings, precautions, and adverse effects.

Contact your Zimmer representative
or visit www.zimmer.com