

Wagner SL Revision[®] Hip Stem

Ordering Information

Implants




Wagner SL Revision
190 mm

H. Wagner

Details
Taper 12/14
Uncemented

Ti6Al7Nb
Protasul[®]-100




Wagner SL Revision
225 mm

H. Wagner

Details
Taper 12/14
Uncemented

Ti6Al7Nb
Protasul-100




Wagner SL Revision
265 mm

H. Wagner

Details
Taper 12/14
Uncemented

Ti6Al7Nb
Protasul-100




Wagner SL Revision
305 mm

H. Wagner

Details
Taper 12/14
Uncemented

Ti6Al7Nb
Protasul-100



STERILE R

Ø mm	REF
14	01.00101.914*
15	01.00101.915*
16	01.00101.916
17	01.00101.917
18	01.00101.918
19	01.00101.919
20	01.00101.920

STERILE R

Ø mm	REF
14	01.00102.214*
15	01.00102.215*
16	01.00102.216
17	01.00102.217
18	01.00102.218
19	01.00102.219
20	01.00102.220
21	01.00102.221
22	01.00102.222

STERILE R

Ø mm	REF
14	01.00102.614*
15	01.00102.615*
16	01.00102.616
17	01.00102.617
18	01.00102.618
19	01.00102.619
20	01.00102.620
21	01.00102.621
22	01.00102.622
23	01.00102.623
24	01.00102.624
25	01.00102.625


STERILE R

Ø mm	REF
14	01.00103.014*
15	01.00103.015*
16	01.00103.016
17	01.00103.017
18	01.00103.018
19	01.00103.019
20	01.00103.020
21	01.00103.021
22	01.00103.022
23	01.00103.023
24	01.00103.024
25	01.00103.025

* Do not use with femoral heads longer than L (+4 mm)

Lengths of 345 mm (34.00.79.XXX) and 385 mm (34.00.79-XXX) available on request and only in the design with CCD of 145°. Impactor/extractor and slot hammer for these selections are also only available by request.

Instruments



Product	REF
Tray basic instruments and awls (complete)	ZS01.00109.000
Tray basic instruments (empty)	01.00109.031
Tray insert basic instruments (empty)	01.00109.032
Tray awls (empty)	01.00109.033
Tray insert awls (empty)	01.00109.034
Standard cover, gray	01.00029.031



Test head for HTP, taper 12/14

Size	Ø mm	REF
S	28	01.01559.128
M	28	01.01559.228
L	28	01.01559.328
XL	28	01.01559.428
S	32	01.01559.132
M	32	01.01559.232
L	32	01.01559.332
XL	32	01.01559.432



Awl for 75.00.25*

Ø mm	REF
14	01.00109.014
15	01.00109.015
16	01.00109.016
17	01.00109.017
18	01.00109.018
19	01.00109.019
20	01.00109.020
21	01.00109.021
22	01.00109.022
23	01.00109.023
24	01.00109.024
25	01.00109.025



Impactor

Product	REF
Handle	75.11.00-02
Synthetic top	75.11.00-03



Handle with quick coupling

REF

75.00.25



Osteotome, straight

Size (mm)	REF
4 × 10	75.11.45-10
5 × 15	75.11.45-15
5 × 20	75.11.45-20



Ruler, long

mm	REF
300	75.11.30



Extractor connector for slap hammer REF 01.00129.190

REF


01.00109.801



Impactor and positioning bar

Product	REF
Impactor	01.00109.808
Positioning bar	75.85.00

* The awl (01.00109.014–025) must under no circumstances be mixed with the awl (75.11.00-140–205).
All instruments

**Product****REF**

Tray trial stems (complete)	ZS01.00109.100
Tray trial stems (empty)	01.00109.035
Standard cover, gray	01.00029.031



Trial stems, distal

Ø mm	REF
14	01.00109.114
15	01.00109.115
16	01.00109.116
17	01.00109.117
18	01.00109.118
19	01.00109.119
20	01.00109.120
21	01.00109.121
22	01.00109.122
23	01.00109.123
24	01.00109.124
25	01.00109.125



Screws for proximal part

mm	REF
190	01.00109.809
225	01.00109.805
265	01.00109.806
305	01.00109.807



Hex wrench

mm	REF
3,5	79.15.84



Trial stems, proximal

mm	REF
190	01.00109.810
225	01.00109.802
265	01.00109.803
305	01.00109.804

Instruments on request



Cerclage pliers

Ø mm	REF
1-2	75.11.00-080



Awls Ø 13 for 75.00.25

REF
01.00109.013



Slot hammer

REF
75.11.00-10



Impactor/extractor
For 345 and 385 mm stems

Ø mm	REF
Basic element	75.11.00-11
Adapter sleeve	75.11.00-12



Test head for THP, taper 12/14

Size	Ø mm	REF
S	36	01.01559.136
M	36	01.01559.236
L	36	01.01559.336
XL	36	01.01559.436



AO Adaptor/TPP Adaptor

REF
5637

X-ray Templates

Description	REF
Wagner SL Revision Stem / 115%	06.01165.000x



Repositioning Tops

Ø mm	REF
28 mm	78.00.38-28
32 mm	78.00.38-32
36 mm	78.00.38-36

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

This material is intended for health care professionals. Distribution to any other recipient is prohibited.

For indications, contraindications, warnings, precautions, potential adverse effects and patient counselling information, see the package insert or contact your local representative; visit www.zimmerbiomet.com for additional product information.

Not for distribution in France. Check for country product clearances and reference product specific instructions for use.

© 2019 Zimmer Biomet

For the Surgical Technique, please refer to 2810.X-GLBL-en Wagner SL Revision Surgical Technique.



0467.1-GLBL-en REV1119



Legal Manufacturer
Zimmer GmbH
Sulzerallee 88404
Winterthur, Switzerland
Telephone +41/ (0)52 262 60 70
Fax +41/ (0)52 262 01 39

www.zimmerbiomet.com



CE mark is not valid unless there is a CE mark on the product label.