



**Acetabular Reconstruction: Cemented Dual Mobility with Acetabular Cups and Cages**

 Indicates "Combination has been assessed as functionally compatible by Zimmer Biomet"  
 Indicates "Combination has either not been assessed or is not functionally compatible"

		Avantage Cemented Cup Size (mm)								
		44	46	48	50	52	54	56	58	60
TMARS Shell Cup Size (mm)	48									
	50									
	52									
	54									
	56									
	58									
	60	✓								
	62	✓	✓							
	64	✓	✓							
	66	✓	✓	✓						
	68	✓	✓	✓						
	70	✓	✓	✓						
	72	✓	✓	✓	✓	✓				
74	✓	✓	✓	✓	✓					
76		✓	✓	✓	✓	✓	✓			
78		✓	✓	✓	✓	✓	✓			
80			✓	✓	✓	✓	✓			

		Avantage Cemented Cup Size (mm)								
		44	46	48	50	52	54	56	58	60
Burch-Scheider Cage Size (mm)	40									
	44									
	50	✓	✓	✓						
	56				✓	✓	✓			
	62							✓	✓	✓

This represents functional compatibility and does not necessarily indicate regulatory approval status for marketing of product pairs.