

- Indicates "Combination has been assessed as functionally compatible by Zimmer Biomet"
- Indicates "Combination has either not been assessed or is not functionally compatible"

Femoral Stems

	G7 Liners Non-Constrained	G7 Liners Constrained
Alloclassic SL and SLO	✓	
Alloclassic SLL Revision Stem	✓	
Alloclassic Variall SLV Stem	✓	
Avenir Cemented	✓	
Avenir Muller Stem	✓	
Avenir Complete Stem	✓	
CLS Sportorno	✓	
CPT 12/14	✓	
CPT 12/14 Long Stem	✓	
CPT 6 degree	✓	
Fitmore Stem	✓	
M/L Taper Stem - Standard and Reduced Neck Length	✓	
M/L Taper with Kinectiv Technology	✓	
Metabloc Cemented	✓	
Metabloc Cementless	✓	
MS-30 stem	✓	
Original M.E. Muller Stem	✓	
TM Primary Hip Stem	✓	
VerSys Advocate Stem	✓	
VerSys CEM/REV/CALCAR (CRC)	✓	
VerSys Fiber Metal Midcoat	✓	
VerSys Fiber Metal Taper	✓	
Versys Beaded Fullcoat (6in 7.5in and 8in)	✓	
Versys Heritage stem	✓	
VerSys LD/Fx cemented	✓	
VerSys LD/Fx Press-Fit	✓	
Wagner Cone Prosthesis Stem	✓	
Wagner SL Revision stem	✓	
CLS Brevious Stem with Kinectiv Technology	✓	
Exafit Stem	✓	
Exafit Revision Stem	✓	
PF Primary Stem	✓	
PF Revision Stem	✓	
Revitan Revision System	✓	
Weber Fix Cemented Stem	✓	
ZMR Revision System	✓	

CoCr Femoral Heads

		G7 Liners Non-Constrained	G7 Liners Constrained
00-8018-028-01	12/14 COCR FEMORAL HEAD 28MM -3.5	✓	
00-8018-028-02	12/14 COCR FEMORAL HEAD 28MM +0	✓	
00-8018-028-03	12/14 COCR FEMORAL HEAD 28MM +3.5	✓	
00-8018-028-14	12/14 COCR FEMORAL HEAD 28MM +7	✓	
00-8018-028-05	12/14 COCR FEMORAL HEAD 28MM +10.5	✓	
00-8018-032-01	12/14 COCR FEMORAL HEAD 32MM -3.5	✓	
00-8018-032-02	12/14 COCR FEMORAL HEAD 32MM +0	✓	
00-8018-032-03	12/14 COCR FEMORAL HEAD 32MM +3.5	✓	
00-8018-032-14	12/14 COCR FEMORAL HEAD 32MM +7	✓	
00-8018-032-05	12/14 COCR FEMORAL HEAD 32MM +10.5	✓	
00-8018-036-01	12/14 COCR FEMORAL HEAD 36MM -3.5	✓	
00-8018-036-02	12/14 COCR FEMORAL HEAD 36MM +0	✓	
00-8018-036-03	12/14 COCR FEMORAL HEAD 36MM +3.5	✓	
00-8018-036-04	12/14 COCR FEMORAL HEAD 36MM +7	✓	
00-8018-036-05	12/14 COCR FEMORAL HEAD 36MM +10.5	✓	
00-8018-040-01	12/14 COCR FEMORAL HEAD 40MM -3.5	✓	
00-8018-040-02	12/14 COCR FEMORAL HEAD 40MM +0	✓	
00-8018-040-03	12/14 COCR FEMORAL HEAD 40MM +3.5	✓	
00-8018-040-04	12/14 COCR FEMORAL HEAD 40MM +7	✓	
00-8018-040-05	12/14 COCR FEMORAL HEAD 40MM +10.5	✓	
00-9026-028-00	6 DEGREE COCR FEMORAL HEAD 28MM +0 (S)	✓	
00-9026-028-35	6 DEGREE COCR FEMORAL HEAD 28MM +3.5 (S+)	✓	
00-9026-029-00	6 DEGREE COCR FEMORAL HEAD 28MM +7 (M)	✓	
00-9026-029-35	6 DEGREE COCR FEMORAL HEAD 28MM +10.5 (M+)	✓	
00-9026-030-00	6 DEGREE COCR FEMORAL HEAD 28MM +14 (L)	✓	
00-9026-032-00	6 DEGREE COCR FEMORAL HEAD 32MM +0 (S)	✓	
00-9026-032-35	6 DEGREE COCR FEMORAL HEAD 32MM +3.5 (S+)	✓	
00-9026-033-00	6 DEGREE COCR FEMORAL HEAD 32MM +7 (M)	✓	
00-9026-033-35	6 DEGREE COCR FEMORAL HEAD 32MM +10.5 (M+)	✓	
00-9026-034-00	6 DEGREE COCR FEMORAL HEAD 32MM +14 (L)	✓	
00-9026-036-00	6 DEGREE COCR FEMORAL HEAD 36MM +0 (S)	✓	
00-9026-036-03	6 DEGREE COCR FEMORAL HEAD 36MM +3.5 (S+)	✓	
00-9026-036-07	6 DEGREE COCR FEMORAL HEAD 36MM +7 (M)	✓	
00-9026-036-10	6 DEGREE COCR FEMORAL HEAD 36MM +10.5 (M+)	✓	
00-9026-036-14	6 DEGREE COCR FEMORAL HEAD 36MM +14 (L)	✓	
01.01002.285	8/10 COCR FEMORAL HEAD 28MM -4 DURASUL	✓	
01.01002.286	8/10 COCR FEMORAL HEAD 28MM +0 DURASUL	✓	
01.01002.287	8/10 COCR FEMORAL HEAD 28MM +4 DURASUL	✓	
01.01012.288	12/14 COCR FEMORAL HEAD 28MM +8 DURASUL	✓	
01.01012.328	12/14 COCR FEMORAL HEAD 32MM +8 DURASUL	✓	
01.01012.365	12/14 COCR FEMORAL HEAD 36MM -4 DURASUL	✓	
01.01012.366	12/14 COCR FEMORAL HEAD 36MM +0 DURASUL	✓	
01.01012.367	12/14 COCR FEMORAL HEAD 36MM +4 DURASUL	✓	
01.01012.368	12/14 COCR FEMORAL HEAD 36MM +8 DURASUL	✓	
01.01012.444	12/14 COCR FEMORAL HEAD 44MM -8mm DURASUL	✓	
01.01012.445	12/14 COCR FEMORAL HEAD 44MM -4mm DURASUL	✓	
01.01012.446	12/14 COCR FEMORAL HEAD 44MM +0 DURASUL	✓	
01.01012.447	12/14 COCR FEMORAL HEAD 44MM +4mm DURASUL	✓	
01.01012.448	12/14 COCR FEMORAL HEAD 44MM +8mm DURASUL	✓	

CoCr Femoral Heads (Cont.)

		G7 Liners Non-Constrained	G7 Liners Constrained

14.28.05-20	12/14 COCR FEMORAL HEAD 28MM -4 (S) DURASUL	✓	
14.28.06-20	12/14 COCR FEMORAL HEAD 28MM +0 (M) DURASUL	✓	
14.28.07-20	12/14 COCR FEMORAL HEAD 28MM +4 (L) DURASUL	✓	
14.32.05-20	12/14 COCR FEMORAL HEAD 32MM -4 (S) DURASUL	✓	
14.32.06-20	12/14 COCR FEMORAL HEAD 32MM +0 (M) DURASUL	✓	
14.32.07-20	12/14 COCR FEMORAL HEAD 32MM +4 (L) DURASUL	✓	
7210-28-000	12/14 COCR FEMORAL HEAD 28MM +0 DURASUL	✓	
7210-28-004	12/14 COCR FEMORAL HEAD 28MM -4 DURASUL	✓	
7210-28-400	12/14 COCR FEMORAL HEAD 28MM +4 DURASUL	✓	
7210-28-800	12/14 COCR FEMORAL HEAD 28MM +8 DURASUL	✓	
7210-32-000	12/14 COCR FEMORAL HEAD 32MM +0 DURASUL	✓	
7210-32-004	12/14 COCR FEMORAL HEAD 32MM -4 DURASUL	✓	
7210-32-400	12/14 COCR FEMORAL HEAD 32MM +4 DURASUL	✓	
7210-32-800	12/14 COCR FEMORAL HEAD 32MM +8 DURASUL	✓	

		G7 Liners Non-Constrained	G7 Liners Constrained
		✓	
Ceramic Femoral Heads			
00-8775-028-01	12/14 BIOLOX DELTA FEMORAL HEAD 28MM -3.5MM	✓	
00-8775-028-02	12/14 BIOLOX DELTA FEMORAL HEAD 28MM +0	✓	
00-8775-028-03	12/14 BIOLOX DELTA FEMORAL HEAD 28MM +3.5	✓	
00-8775-032-01	12/14 BIOLOX DELTA FEMORAL HEAD 32MM -3.5	✓	
00-8775-032-02	12/14 BIOLOX DELTA FEMORAL HEAD 32MM +0	✓	
00-8775-032-03	12/14 BIOLOX DELTA FEMORAL HEAD 32MM +3.5	✓	
00-8775-032-04	12/14 BIOLOX DELTA FEMORAL HEAD 32MM +7	✓	
00-8775-036-01	12/14 BIOLOX DELTA FEMORAL HEAD 36MM -3.5	✓	
00-8775-036-02	12/14 BIOLOX DELTA FEMORAL HEAD 36MM +0	✓	
00-8775-036-03	12/14 BIOLOX DELTA FEMORAL HEAD 36MM +3.5	✓	
00-8775-036-04	12/14 BIOLOX DELTA FEMORAL HEAD 36MM +7	✓	
00-8775-040-01	12/14 BIOLOX DELTA FEMORAL HEAD 40MM -3.5	✓	
00-8775-040-02	12/14 BIOLOX DELTA FEMORAL HEAD 40MM +0	✓	
00-8775-040-03	12/14 BIOLOX DELTA FEMORAL HEAD 40MM +3.5	✓	
00-8775-040-04	12/14 BIOLOX DELTA FEMORAL HEAD 40MM +7	✓	
00-8777-028-01	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 28MM -3.0	✓	
00-8777-028-02	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 28MM +0	✓	
00-8777-028-03	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 28MM +3.5	✓	
00-8777-028-04	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 28MM +7	✓	
00-8777-032-01	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 32MM -3.0	✓	
00-8777-032-02	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 32MM +0	✓	
00-8777-032-03	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 32MM +3.5	✓	
00-8777-032-04	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 32MM +7	✓	
00-8777-036-01	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 36MM -3.0	✓	
00-8777-036-02	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 36MM +0	✓	
00-8777-036-03	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 36MM +3.5	✓	
00-8777-036-04	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 36MM +7	✓	
00-8777-040-01	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 40MM -3.0	✓	
00-8777-040-02	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 40MM +0	✓	
00-8777-040-03	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 40MM +3.5	✓	
00-8777-040-04	12/14 BIOLOX DELTA OPTION FEMORAL HEAD 40MM +7	✓	
00-8778-028-01	8/10 BIOLOX OPTION Hd/Adpt 28MM -3.0	✓	
00-8778-028-02	8/10 BIOLOX OPTION Hd/Adpt 28MM +0	✓	
00-8778-028-03	8/10 BIOLOX OPTION Hd/Adpt 28MM +4.0	✓	
00-8778-032-01	8/10 BIOLOX OPTION Hd/Adpt 32MM -3.0	✓	
00-8778-032-02	8/10 BIOLOX OPTION Hd/Adpt 32MM +0	✓	
00-8778-032-03	8/10 BIOLOX OPTION Hd/Adpt 32MM +4.0	✓	
00-8778-036-01	8/10 BIOLOX OPTION Hd/Adpt 36MM -3.0	✓	
00-8778-036-02	8/10 BIOLOX OPTION Hd/Adpt 36MM +0	✓	
00-8778-036-03	8/10 BIOLOX OPTION Hd/Adpt 36MM +4.0	✓	
00-8778-040-01	8/10 BIOLOX OPTION Hd/Adpt 40MM -3.0	✓	
00-8778-040-02	8/10 BIOLOX OPTION Hd/Adpt 40MM +0	✓	
00-8778-040-03	8/10 BIOLOX OPTION Hd/Adpt 40MM +4.0	✓	
00-6428-028-01	12/14 ALUMINA CERAMIC FEMORAL HEAD 28MM -3.5	✓	
00-6428-028-02	12/14 ALUMINA CERAMIC FEMORAL HEAD 28MM +0	✓	
00-6428-028-03	12/14 ALUMINA CERAMIC FEMORAL HEAD 28MM +3.5	✓	
00-6428-032-01	12/14 ALUMINA CERAMIC FEMORAL HEAD 32MM -3.5	✓	
00-6428-032-02	12/14 ALUMINA CERAMIC FEMORAL HEAD 32MM +0	✓	
00-6428-032-03	12/14 ALUMINA CERAMIC FEMORAL HEAD 32MM +3.5	✓	
12.28.05	12/14 BIOLOX FORTE FEMORAL HEAD 28MM -3.5	✓	
12.28.06	12/14 BIOLOX FORTE FEMORAL HEAD 28MM +0	✓	
12.28.07	12/14 BIOLOX FORTE FEMORAL HEAD 28MM +3.5	✓	
12.32.05	12/14 BIOLOX FORTE FEMORAL HEAD 32MM -4	✓	
12.32.06	12/14 BIOLOX FORTE FEMORAL HEAD 32MM +0	✓	
12.32.07	12/14 BIOLOX FORTE FEMORAL HEAD 32MM +4	✓	
Stainless Steel Femoral Heads			
		G7 Liners Non-Constrained	G7 Liners Constrained
		✓	

00-8011-001-28	6 DEGREE CPT SST FEMORAL HEAD 28MM SHT NECK	✓	
00-8011-001-32	6 DEGREE CPT SST FEMORAL HEAD 32MM SHT NECK	✓	
00-8011-002-28	6 DEGREE CPT SST FEMORAL HEAD 28MM MED NECK	✓	
00-8011-002-32	6 DEGREE CPT SST FEMORAL HEAD 32MM MED NECK	✓	
00-8011-003-28	6 DEGREE CPT SST FEMORAL HEAD 28MM LNG NECK	✓	
00-8011-003-32	6 DEGREE CPT SST FEMORAL HEAD 32MM LNG NECK	✓	
01.01102.285	8/10 PROTASUL S-30 FEMORAL HEAD 28MM -4 (S)	✓	
01.01102.286	8/10 PROTASUL S-30 FEMORAL HEAD 28MM +0 (M)	✓	
01.01102.287	8/10 PROTASUL S-30 FEMORAL HEAD 28MM +4 (L)	✓	
30.28.05	12/14 PROTASUL S-30 FEMORAL HEAD 28MM -4 (S)	✓	
30.28.06	12/14 PROTASUL S-30 FEMORAL HEAD 28MM +0 (M)	✓	
30.28.07	12/14 PROTASUL S-30 FEMORAL HEAD 28MM +4 (L)	✓	
30.32.05	12/14 PROTASUL S-30 FEMORAL HEAD 32MM -4 (S)	✓	
30.32.06	12/14 PROTASUL S-30 FEMORAL HEAD 32MM +0 (M)	✓	
30.32.07	12/14 PROTASUL S-30 FEMORAL HEAD 32MM +4 (L)	✓	