

**Articulating Combinations:  
Ceramic on Ceramic**

- Indicates "Combination is approved by Zimmer, Inc."
- Indicates "Combination is not approved by Zimmer, Inc."

**Ceramic Femoral Heads:**

- Alumina 12/14 Femoral Head (XX-6418/6428-XXX-XX)
- BIOLOX-forte head (12.28.32.05/07)
- CERASUL head
- Continuum Alumina Bearing Head (XX-4102-XXX-XX)
- Continuum Zirconia Bearing Head (XX-4101-XXX-XX)
- ESKA BIOLOX Forte head
- NGK Zirconia femoral head
- SULOX head
- Tete Alumine Cone 8/10 D28
- Zirconia head (7676-XX-XXX)
- BIOLOX delta head (00-8775-028\_040-XX)(12/14 taper)
- BIOLOX OPTION head system (00-8777-028\_040-XX)(12/14 taper)
- BIOLOX OPTION head system (00-8778-028\_36-XX)(8/10 taper)
- BIOLOX OPTION head system (00-8778-044\_048-XX) (12/14 taper)

**Ceramic Inserts (Liners):**

- ALPHA CERASUL insert
- BIOLOX Delta Taper Liner (8775)\*\*
- ESKA Spongiosa alumina insert
- GAMMA CERASUL insert
- Trilogy AB Alumina Insert \*
- ALLOCLASSIC Variall Biolox Delta Insert

		✓								✓	✓	✓		
					✓									
		✓								✓	✓	✓		
✓										✓	✓	✓		
										✓	✓	✓		

**Ceramic on Ceramic Monoblock Cups:**

- Maxera Cup

									✓	✓	✓	✓
--	--	--	--	--	--	--	--	--	---	---	---	---

\* Trilogy AB Alumina Insert with BIOLOX delta and OPTION heads: This pairing is not released in the U.S.!

\*\* Limited to investigational use in the United States.