

Zimmer Intertrochanteric/Subtrochanteric Fixation Devices (ITST) System Compatibility

- Indicates that the combination CAN be used
- Indicates that the combination CANNOT be used

System Designation	ITST	M/DN	ITST	ITST	ITST	ITST
<p>Note: This table represents functional compatibility and does not indicate the regulatory approval status.</p> <p>Be sure to check for proper Regulatory approvals in your country prior to using any products found on this website.</p> <p>Some devices may not be currently licensed with Health Canada. Some device compatibilities may not be approved for use by Health Canada.</p> <p>xx or xxx = These are placeholders that represent the item numbers for the range of device sizes.</p>	ITST 4.5mm Diameter Cortical Screw, 20-90mm Length 00-2253-0xx-45	M/DN 5.5mm Diameter Cortical Screw, 20-100mm Length 00-2253-xxx-55	ITST One-Piece Lag Screw Asia, 11mm Diameter, 70-120mm Length 00-2256-002-xx	ITST One-Piece Lag Screw, 11mm Diameter, 70-120mm Length 00-2259-001-xx	ITST Anti-Rotation Screw, 6.5mm, 60-110mm Length 00-2257-xxx-65	ITST One Piece Nail Cap, 0-20mm Head Height 00-2259-007/008/009-xx
ITST Femoral Nail, 10mm Diameter, 18cm Length, Asia (Left/Right) 00-2256-180-10	✓	✗	✓	✗	✓	✓
ITST Femoral Nail, 11-15mm Diameter, 18cm Length, Asia (Left/Right) 00-2256-180-11/12/13/14/15	✗	✓	✓	✗	✓	✓
ITST Femoral Nail, 10mm Diameter, 18-50cm Length (Left/Right) 00-2257-xx0-00/10	✓	✗	✗	✓	✓	✓
ITST Femoral Nail, 11-15mm Diameter, 18-50cm Length (Left/Right) 00-2257-xx0-11/12/13/14/15 00-2257-xx0-01/02/03/04	✗	✓	✗	✓	✓	✓