

# Zimmer M/DN System Compatibility

Revised 9/22/2015



Indicates that the combination CAN be used

Indicates that the combination CANNOT be used

System Designation	M/DN	M/DN	M/DN	M/DN	M/DN	M/DN	M/DN
<p>Note: This table represents functional compatibility and does not indicate the regulatory approval status.</p> <p><b>Be sure to check for proper Regulatory approvals in your country prior to using any products found on this website.</b></p> <p>Some devices may not be currently licensed with Health Canada. Some device compatibilities may not be approved for use by Health Canada.</p> <p>xx or xxx = These are placeholders that represent the item numbers for the range of device sizes.</p>	M/DN 3.7mm Diameter Cortical Screw, 20-90mm Length 00-2253-0xx-37	M/DN 4.2mm Diameter Cortical Screw, 20-90mm Length 00-2253-0xx-42	M/DN 5.5mm Diameter Cortical Screw, 20-100mm Length 00-2253-xxx-55	M/DN 5.5mm Dia., 3.5mm Hex Head Recon Screw, 60-130mm Length 00-2252-xxx-55	M/DN Extended Nail End Cap, 0-15mm Height 00-2252-001-xx	M/DN Humeral Locking Nail Cap 00-2252-003-00	M/DN Retrograde Nail Locking Cap 00-2252-004-00
M/DN Retrograde Femoral IM Nail, 9-11mm Diameter, 14-48cm Length 00-2240-009/010/011-xx		✓	✓				✓
M/DN Retrograde Femoral IM Nail, 12-14mm Diameter, 14-48cm Length 00-2240-012/013/014-xx			✓				✓
M/DN Retrograde 5.5mm Diameter Femoral Cortical Lock Nut 00-2240-050-00			✓				
M/DN Femoral Recon Nail, 8mm Diameter, 24-44cm 00-2252-xx0-08	✓		✓	✓	✓		
M/DN Femoral Recon Nail, 9-11mm Diameter, 24-48cm 00-2252-xx0-09/010/011		✓	✓	✓	✓		
M/DN Femoral Recon Nail, 12-16mm Diameter, 24-50cm 00-2252-xx0-012/013/014/015/016			✓	✓	✓		
M/DN 13mm Washer 00-2252-002-00			✓				
M/DN Tibial Nail, 6mm, 9-15mm Diameter, 18-46mm Length 00-2253-xx0-06/09/10/11/12/13/14/15		✓			✓		
M/DN Tibial Nail, 7-8mm Diameter, 18-42mm Length 00-2253-xx0-07/08	✓	✓			✓		
M/DN Humeral Nail, 6mm, 9-13mm Dia., 18-30mm Length 00-2254-xxx-06/09/10/11/12/13		✓				✓	
M/DN Humeral Nail, 7-8mm Diameter, 18-30mm Length 00-2254-xxx-07/08	✓	✓				✓	

Note: This table represents functional compatibility and does not indicate the regulatory approval status.