
 Indicates that the combination CAN be used
 Indicates that the combination CANNOT be used

Screw Head Size System Designation		2.7mm			4.0mm			3.5mm					4.5mm										
		ZPS	ZPS	Peri	Peri	Peri	Peri	ZPS	ZPS	ZPS	ZPS	ZPS	Peri	Peri	Peri	Peri	ZPS	ZPS	ZPS	ZPS	ZPS	ZPS	
Note: This table represents functional compatibility and does not indicate the regulatory approval status. xx or xxx = These are placeholders that represent the item numbers for the range of device sizes.		2.7mm Cortical Screw, Self-Tapping 47/00-48-27-xxx-01 2.7mm Cortical Screw 47/00-48-27-xxx-00 3.5mm Periarticular Cortical Screw with 2.7mm Head, Self-Tapping 47/00-23-48-xxx-35 4.0mm Forte Cancellous Screw, Partially Threaded 47/00-24-80-xxx-40 4.0mm Forte Cancellous Screw, Fully Threaded 47/00-24-80-xxx-41 3.5mm Periarticular Cortical Screw, Self-Tapping 47/00-23-47-020/120-xx 4.0mm Periarticular Cancellous Screw, Partially Threaded 47/00-23-47-021-xx 4.0mm Periarticular Cancellous Screw, Fully Threaded 47/00-23-47-022-xx 3.5mm Cortical Screw 47/00-48-35-xxx-00 3.5mm Cortical Screw, Self-Tapping 47/00-48-35-xxx-01 3.5mm Cortical Shaft Screw 47/00-48-35-xxx-04 4.0mm Cancellous Screw, Fully Threaded 47/00-48-60-xxx-00 4.0mm Cancellous Screw, Partially Threaded 47/00-48-60-xxx-01 4.5mm Periarticular Cortical Screw, Self-Tapping 47/00-23-47-023/123-xx 6.5mm Periarticular Cancellous Screw, 16mm Thread 47/00-23-47-024/124-xx 6.5mm Periarticular Cancellous Screw, 32mm Thread 47/00-23-47-025/125-xx 4.5mm Periarticular Cancellous Screw, Fully Threaded 47/00-23-47-026/126-xx 4.5mm Cortical Screw 47/00-48-45-xxx-00 4.5mm Cortical Screw, Self-Tapping 47/00-48-45-xxx-01 4.5mm Cortical Shaft Screw 47/00-48-45-xxx-04 4.5mm Multiaxial Screw 47/00-48-45-xxx-08 6.5mm Cancellous Screw, Fully Threaded 47/00-48-65-xxx-00 6.5mm Cancellous Screw, 16mm Thread 47/00-48-65-xxx-01 6.5mm Cancellous Screw, 32mm Thread 47/00-48-65-xxx-02																					
3.5mm Plates	Periarticular Plating System 3.5mm Proximal Lateral Humeral Plate (Left/Right) 00-2348-001/002-xx																						
	Periarticular Plating System 3.5mm Proximal Tibial Plates (Left/Right) 00-2347-003/004/011/012-xx																						
4.5mm Plates	Periarticular Plating System 4.5mm Proximal Lateral Humeral Plate (Left/Right) 00-2348-003/004-xx																						
	Periarticular Plating System 4.5mm Distal Medial Femoral Plate (Left/Right) 00-2347-031/032-xx																						
	Periarticular Plating System 4.5mm Proximal Tibial Plate (Left/Right) 00-2347-005/006/013/014-xx																						
Humeral Plates	Periarticular Plating System Distal Humeral Plates (Left/Right) 00-2348-005/006/007-xx																						
Ulnar Plates	Periarticular Plating System Proximal Dorsal Ulnar Plate (Left/Right) 00-2348-009/010-xx																						
Radial Plates	Periarticular Plating System Small Distal Radial Plate (Left/Right) 00-2348-013/014/017/018/021/022/025/026-xx																						
Femoral Plates	Periarticular Plating System Distal Lateral Femoral Plate (Left/Right) 00-2347-001/002-xx																						
Tibial Plates	Periarticular Plating System Distal Lateral/Medial Tibial Plate (Left/Right) 00-2347-007/008/009/010-xx																						
Fibular Plates	Periarticular Plating System Distal Posterior/Lateral Fibular Plate (Left/Right) 00-2347-015/016/017/018-xx																						
Calcaneal Plates	Periarticular Plating System Lateral Calcaneal Plate 00-2347-035/036-00																						

Note: This table represents functional compatibility and does not indicate the regulatory approval status.