



**Zimmer Versa-Fx/Versa-Fx II Femoral Fixation System and
Free-Lock Femoral Hip Fixation System Compatibility**

 Indicates that the combination is approved by Zimmer
 Indicates that the combination IS NOT approved by Zimmer

System Designation		Degree Angle	Screw Head Size		System Designation		Lag Screws		Compression Screws			
			4.5mm heads	6.5mm heads	ZPS	ZPS	Magna-FX	Versa-Fx	Versa-Fx			
Versa-FX II		90° Angle	<p>Note: This table represents functional compatibility and does not indicate the regulatory approval status.</p> <p>Be sure to check for proper Regulatory approvals in your country prior to using any products found on this website.</p> <p>Some devices may not be currently licensed with Health Canada. Some device compatibilities may not be approved for use by Health Canada.</p> <p>xx or xxx = These are placeholders that represent the item numbers for the range of device sizes.</p>									
			4.5mm Cortical Screws 47/00-4845-xxx-00	4.5mm Cortical Screws, Self-Taping 47/00-4845-xxx-01	6.5mm Cancellous Screws (16mm Threaded) 47/00-4865-xxx-01	6.5mm Cancellous Screws (32mm Threaded) 47/00-4865-xxx-02	6.5mm Cancellous Screws (Fully Threaded) 47/00-4865-xxx-00	7.0mm Cannulated Bone Screws (16mm Threaded) 00-1146-xxx-00	7.0mm Cannulated Bone Screws (32mm Threaded) 00-1146-xxx-32	7.0mm Cannulated Bone Screws (Fully Threaded) 00-1146-xxx-99	1193 Series Versa-Fx Lag Screws (Short) 00-1193-005-xx	1193 Series Versa-Fx Lag Screws (Standard) 00-1193-010-xx
Versa-FX		90° Angle	Supracondylar Tube/Plate Low Profile 00-1193-090-xx	✓	✓	✓	✓	✓	✓	✓	✓	✓
		95° Angle	Supracondylar Tube/Plate Low Profile 00-1193-095-xx	✓	✓	✓	✓	✓	✓	✓	✓	✓
Versa-FX II		130° Angle	Proximal Femur Keyed Compression Tube/Plate Standard Low Profile 00-1193-130-xx	✓	✓	✓	✓	✓	✓	✓	✓	✓
			Proximal Femur Keyed Compression Tube/Plate Short Low Profile 00-1193-130-94/95/96	✓	✓	✓	✓	✓	✓	✓	✓	✓
Versa-FX II		135° Angle	Proximal Femur Compression Tube/Plate Standard Low Profile - Keyless 00-1194-130-04/05/06/08	✓	✓	✓	✓	✓	✓	✓	✓	✓
			Proximal Femur Keyed Compression Tube/Plate Standard Low Profile 00-1193-135-xx	✓	✓	✓	✓	✓	✓	✓	✓	✓
Versa-FX II		140° Angle	Proximal Femur Keyed Compression Tube/Plate Short Low Profile 00-1193-135-94/95/96	✓	✓	✓	✓	✓	✓	✓	✓	✓
			Proximal Femur Compression Tube/Plate Standard Low Profile - Keyless 00-1194-135-04/05/06/08	✓	✓	✓	✓	✓	✓	✓	✓	✓
Versa-FX II		145° Angle	Proximal Femur Keyed Compression Tube/Plate Standard Low Profile 00-1193-140-xx	✓	✓	✓	✓	✓	✓	✓	✓	✓
			Proximal Femur Keyed Compression Tube/Plate Short Low Profile 00-1193-140-94/95/96	✓	✓	✓	✓	✓	✓	✓	✓	✓
Versa-FX II		150° Angle	Proximal Femur Compression Tube/Plate Standard Low Profile - Keyless 00-1194-140-04/05/06/08	✓	✓	✓	✓	✓	✓	✓	✓	✓
			Proximal Femur Keyed Compression Tube/Plate Standard Low Profile 00-1193-145-xx	✓	✓	✓	✓	✓	✓	✓	✓	✓
Versa-FX		90° Angle	Proximal Femur Keyed Compression Tube/Plate Short Low Profile 00-1193-145-94/95/96	✓	✓	✓	✓	✓	✓	✓	✓	✓
			Proximal Femur Compression Tube/Plate Standard Low Profile - Keyless 00-1194-145-04/05/06/08	✓	✓	✓	✓	✓	✓	✓	✓	✓
Free-Lock		130° Angle	Proximal Femur Keyed Compression Tube/Plate Standard Low Profile 00-1193-145-xx	✓	✓	✓	✓	✓	✓	✓	✓	✓
			Proximal Femur Keyed Compression Tube/Plate Short Low Profile 00-1193-145-94/95/96	✓	✓	✓	✓	✓	✓	✓	✓	✓
Free-Lock		135° Angle	Proximal Femur Keyed Compression Tube/Plate Standard Low Profile 00-1193-150-xx	✓	✓	✓	✓	✓	✓	✓	✓	✓
			Proximal Femur Keyed Compression Tube/Plate Short Low Profile 00-1193-150-94/95/96	✓	✓	✓	✓	✓	✓	✓	✓	✓
Free-Lock		140° Angle	Proximal Femur Keyed Compression Tube/Plate Standard Low Profile 00-1193-150-xx	✓	✓	✓	✓	✓	✓	✓	✓	✓
			Proximal Femur Keyed Compression Tube/Plate Short Low Profile 00-1193-150-94/95/96	✓	✓	✓	✓	✓	✓	✓	✓	✓
Free-Lock		145° Angle	Proximal Femur Keyed Compression Tube/Plate Standard Low Profile 00-1193-150-xx	✓	✓	✓	✓	✓	✓	✓	✓	✓
			Proximal Femur Keyed Compression Tube/Plate Short Low Profile 00-1193-150-94/95/96	✓	✓	✓	✓	✓	✓	✓	✓	✓
Free-Lock		150° Angle	Proximal Femur Keyed Compression Tube/Plate Standard Low Profile 00-1193-150-xx	✓	✓	✓	✓	✓	✓	✓	✓	✓
			Proximal Femur Keyed Compression Tube/Plate Short Low Profile 00-1193-150-94/95/96	✓	✓	✓	✓	✓	✓	✓	✓	✓

Compression Screws to Lag Screw Compatibility	
System Designation	Compatibility
Versa-Fx Lag Screw	1193 Series Versa-Fx Lag Screws (Short) 00-1193-005-xx
Versa-Fx Lag Screw	1193 Series Versa-Fx Lag Screws (Standard) 00-1193-010-xx
Versa-Fx Lag Screw	1193 Series Versa-Fx Lag Screws (Large) 00-1193-015-xx
Versa-Fx Lag Screw	4.7mm Hex Head Compression Screws 47-1162-xxx-00

Note: This table represents functional compatibility and does not indicate the regulatory approval status.