

Zimmer Periarticular Locking Plate (ZPLP) System Compatibility

- ✓ Indicates that the combination CAN be used
- ✗ Indicates that the combination CANNOT be used

Screw Head Size	2.4mm heads			2.7mm heads					3.5mm heads										5.5mm heads										3.5mm	5.5mm									
	System Designation	ZPLP	ZPLP	ZPLP	ZPLP	ZPLP	ULS	Peri	ZPS	ZPLP	ZPLP	ZPLP	ZPLP	ZML	Peri	Peri	Peri	ZPS	ZPS	ZPS	ZPS	ZPS	ZPS	ZPLP	ZPLP	ZPLP	ZML	Peri	Peri	Peri	Peri	ZPS	ZPS	ZPS	ZPS	ZPS	ZPS	Cable-Ready	Cable-Ready
<p>Note: This table represents functional compatibility and does not indicate the regulatory approval status.</p> <p>Be sure to check for proper Regulatory approvals in your country prior to using any products found on this website.</p> <p>Some devices may not be currently licensed with Health Canada. Some device compatibilities may not be approved for use by Health Canada.</p> <p>xx or xxx = These are placeholders that represent the item numbers for the range of device sizes.</p>																																							
<p>1.8mm Locking Peg 47/00-2359-xxx-18</p> <p>2.4mm Locking Screw 47/00-2359-xxx-24</p> <p>2.4mm Conical Screw 47/00-2359-xxx-25</p> <p>2.7mm Conical Screw 47/00-2359-xxx-26</p> <p>3.5mm Locking Screw with 2.7mm Head 47/00-2359-xxx-38</p> <p>2.7mm Locking Screw 47/00-4828-xxx-02</p> <p>3.5mm Periarticular Screw with 2.7mm Head 47/00-2348-xxx-35</p> <p>2.7mm Cortical Screw, Self-Tapping 47/00-4827-xxx-01</p> <p>3.5mm locking Screw 47/00-2359-xxx-35</p> <p>3.5mm Cannulated Locking Screw 47/00-2359-xxx-36</p> <p>3.5mm Cannulated Conical Screw 47/00-2359-xxx-37</p> <p>2.7mm Locking Screw 47/00-2359-xxx-27</p> <p>3.5mm Zimmer MotionLoc Screw 47/00-2359-xxx-39</p> <p>3.5mm Periarticular Screw, Self-Tapping 47/00-2347-020/120-xx</p> <p>4.0mm Periarticular Screw, Partially Threaded 47/00-2347-021-xx</p> <p>4.0mm Periarticular Screw, Fully Threaded 47/00-2347-022-xx</p> <p>3.5mm Cortical Screw 47/00-4835-xxx-00</p> <p>3.5mm Cortical Screw, Self-Tapping 47/00-4835-xxx-01</p> <p>3.5mm Cortical Shaft Screw 47/00-4835-xxx-04</p> <p>4.0mm Cancellous Screw, Fully Threaded 47/00-4840-xxx-00</p> <p>4.0mm Cancellous Screw, Partially Threaded 47/00-4840-xxx-01</p> <p>4.5mm Locking Screw 47/00-2359-xxx-45</p> <p>5.5mm Cannulated Locking Screw 47/00-2359-xxx-55</p> <p>5.5mm Cannulated Conical Screw 47/00-2359-xxx-56</p> <p>4.5mm Zimmer MotionLoc Screw 47/00-2359-xxx-49</p> <p>4.5mm Periarticular Screw, Self-Tapping 47/00-2347-023/123-xx</p> <p>6.5mm Periarticular Cancellous Screw, 16mm Thread 00-2347-024/124-xx</p> <p>6.5mm Periarticular Cancellous Screw, 32mm Thread 00-2347-025/125-xx</p> <p>6.5mm Periarticular Cancellous Screw, Fully Threaded 00-2347-026/126-xx</p> <p>4.5mm Cortical Screw 47/00-4845-xxx-00</p> <p>4.5mm Cortical Screw, Self-Tapping 47/00-4845-xxx-01</p> <p>4.5mm Cortical Shaft Screw 47/00-4845-xxx-04</p> <p>6.5mm Cancellous Screw, Fully Threaded 47/00-4865-xxx-00</p> <p>6.5mm Cancellous Screw, 16mm thread 47/00-4865-xxx-01</p> <p>6.5mm Cancellous Screw, 32mm Thread 47/00-4865-xxx-02</p> <p>Cable-Ready Cable Button for ZPLP 3.5mm Locking Screw Holes 47-2232-060-35</p> <p>Cable-Ready Cable Button for ZPLP 5.5mm Locking Screw Holes 47-2232-060-55</p>																																							
2.4mm plates	Periarticular Locking Plate Distal Radial Locking Plates (Left/Right) 47/00-2358-013/014/017/018/021/022/025/026-xx	✓	✓	✓				✓	✓																														
	Periarticular Locking Plate Radial Styloid Locking Plates (Left/Right) 47/00-2358-027/028-xx	✓	✓	✓					✓																														
2.7mm plates	Periarticular Locking Plate Distal Lateral Fibular Locking Plates (Left/Right) 47/00-2357-017/018-xx				✓	✓	✓	✓	✓																														
	Periarticular Locking Plate Distal Humeral and Proximal Ulnar Locking Plates (Left/Right) 47/00-2358-005/006/007/008/107/108/009/012/023/024-xx				✓	✓	✓	✓	✓																														
3.5mm plates	Periarticular Locking Plate Proximal Lateral Humeral Locking Plates (Left/Right) 47/00-2358-001/002-xx									✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓															✓	
	Periarticular Locking Plate Distal Tibial Locking Plates (Left/Right) 47/00-2357-007/008/009/010-xx									✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓																✓	
	Periarticular Locking Plate Proximal Tibial Locking Plates (Left/Right) 47/00-2357-003/004-xx									✓	✓	✓		✓	✓	✓	✓	✓	✓	✓	✓	✓																✓	
5.5mm plates	Periarticular Locking Plate Proximal Lateral Tibial Locking Plates (Left/Right) 47/00-2357-005/006-xx																						✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Periarticular Locking Plate Distal Lateral Femoral Locking Plates (Left/Right) 47/00-2357-101/102-xx																						✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

This represents functional compatibility and does not necessarily indicate regulatory approval status for marketing of product pairs.